

Telligence 7 Elevates the Clinical Workflow Experience

Simplifying technology ownership with 7 new features



Cloud-based installation

Cut IT installation time, reduce costs and increase cyber security and scalability using a SQL server



Unmatched efficiency with Import and Export templates

Increase installation speed and get service time on site down to just 30 minutes



No downtime for IP restarts

Make changes and apply them in real time without any disruption to your staff



Say hello to the new room gateway

This new hardware saves space, reduces cost and supports expanded capabilities, like medical device integration



Paging just got easier

With Ascom's PBX integration through a SIP connection, simply use the handset on the nurse call console for paging



No more ripping and replacing cable

Spend less to get more when you move to Ascom Telligence nurse call by using your existing cabling



Ultimate customization

Modify all buttons on the bedside modules, from system wide to room level. Integrate connections to virtual nursing platforms with our new pillow speaker model

Secure, adaptable and affordable

For your hospital or Enterprise requirements and deployments